



nano

CUSTOM
FIT





PERFECT FIT. ABSOLUTE COMFORT.
TOTAL ADAPTABILITY.

Ready for the new era of Nano glasses?
Introducing **Nano Custom Fit 3.0!**

Made from Siliflex™—our patented material that makes Nano Vista glasses ultra-resistant, super flexible, lightweight, and hypoallergenic. They always retain their original shape, even when bent or twisted, because they are truly **#indestructible!**

#indestructible

Nano Custom Fit takes it a step further...
adapting like never before!

Designed for kids and teens aged 4 to 16,
featuring Hinge 3.0 and adjustable temple
tips for superior comfort.

Includes both full strap and mini strap.



And as always, thinking of opticians!

All models include demo lenses with a measurement grid
to simplify PD and fitting height measurements.



THE SECRET BEHIND THE PERFECT FIT: INNOVATION THAT ADAPTS TO EVERY CHILD!



Geometry adapted with Artificial Intelligence

Our best-selling Nano Indestructible models, now with redesigned front geometry to fit 100% of children around the world.




Extendable Temples: A Perfect Fit for Every Child

An innovative system that allows adjustment of temple length to suit different ear-to-eye positions, ensuring an optimal fit and maximum stability for every user.



Adjustable Nose Pad KIU4: Guaranteed Comfort

New design with a flexible metal core and a deeper, rounded surface offering multi-directional angle and width adjustment for a precise and comfortable fit across a range of nasal anatomies.



NANO CUSTOM FIT: THE EVOLUTION OF INDESTRUCTIBLE GLASSES!



REPLAY CF

CALIBER



44 (4-6 years) 44mm 112-127mm 13mm 31mm
46 (6-8 years) 46mm 125-140mm 14mm 32,5mm

Measures according to the DIN EN ISO 8624 standard.



BLACK/BLUE

44: **NAO3320144** 46: **NAO3320146**



CRYSTAL PURPLE/PINK

44: **NAO3320344** 46: **NAO3320346**



CRYSTAL SOFT BLUE TURQUOISE/BLUE

44: **NAO3320544** 46: **NAO3320546**



MATTE NAVY/ORANGE

44: **NAO3320244** 46: **NAO3320246**



MATTE SOFT PINK/PINK

44: **NAO3320444** 46: **NAO3320446**





nano

CUSTOM
FIT





nano

CUSTOM
FIT



F L I C K E R_{CF}

CALIBER



44 (4-6 years)	44mm	112-127mm	15mm	37,5mm
46 (6-8 years)	46mm	125-140mm	16mm	39mm

Measures according to the DIN EN ISO 8624 standard.



BLACK/BLACK

44: **NAO3380144** 46: **NAO3380146**



MATTE NAVY/BUE

44: **NAO3380244** 46: **NAO3380246**



CRYSTAL SOFT BLUE/BLUE

44: **NAO3380444** 46: **NAO3380446**



CRYSTAL GARNET/GREY

44: **NAO3380344** 46: **NAO3380346**



CRYSTAL SOFT PINK/PINK

44: **NAO3380544** 46: **NAO3380546**





CALIBER    

46 (6-8 years)	46mm	125-140mm	15mm	31,5mm
48 (8-12 years)	48mm	125-140mm	16mm	33mm

Measures according to the DIN EN ISO 8624 standard.



MATTE PINK/LIGHT PINK

46: **NAO3330246** 48: **NAO3330248**



MATTE BLACK/RED

46: **NAO3330146** 48: **NAO3330148**



CRYSTAL NAVY/GREEN

46: **NAO3330346** 48: **NAO3330348**



CRYSTAL SOFT BLUE/BLUE

46: **NAO3330446** 48: **NAO3330448**



CRYSTAL LILAC/PURPLE

46: **NAO3330546** 48: **NAO3330548**





nano
CUSTOM
FIT



CAMP ER CF

CALIBER    

46 (6-8 years) 46mm 125-140mm 15mm 33,5mm
48 (8-12 years) 48mm 125-140mm 16mm 35mm

Measures according to the DIN EN ISO 8624 standard.



BLACK/FUCHSIA

46: **NAO3350146** 48: **NAO3350148**



CRYSTAL RASPBERRY/RASPBERRY

46: **NAO3350346** 48: **NAO3350348**



CRYSTAL NAVY/BLUE/ORANGE

46: **NAO3350246** 48: **NAO3350248**



MATTE SOFT LILAC/PURPLE

46: **NAO3350446** 48: **NAO3350448**



SOFT BLUE/BLUE

46: **NAO3350546** 48: **NAO3350548**






nano
CUSTOM
FIT



FANGAME_{CF}

CALIBER



50 (12-14 years) 50mm 125-140mm 16mm 34mm
52 (14-16 years) 52mm 125-140mm 16mm 35mm

Measures according to the DIN EN ISO 8624 standard.



MATTE GREY/DARK GREY/RED

50: NAO3340250 52: NAO3340252



CRYSTAL PINK/PINK

50: NAO3340450 52: NAO3340452



BLACK/BLACK

50: NAO3340150 52: NAO3340152



MATTE NAVY/BLUE


50: NAO3340350 52: NAO3340352



CRYSTAL PEARL GREY/GREY/BLUE

50: NAO3340550 52: NAO3340552





nano
CUSTOM
FIT





nano
CUSTOM
FIT



MULTIPLAYER_{CF}



BLACK/BLUE
48: **NAO3390148** 50: **NAO3390150**



MATTE NAVY/GREEN
48: **NAO3390248** 50: **NAO3390250**



CRYSTAL LIGHT GREY/BLUE
48: **NAO3390448** 50: **NAO3390450**



MATTE GREY/RED/BLUE
48: **NAO3390348** 50: **NAO3390350**



CRYSTAL PURPLE/DARK PURPLE
48: **NAO3390548** 50: **NAO3390550**

CALIBER



48 (8-12 years)	48mm	125-140mm	18mm	43mm
50 (12-14 years)	50mm	125-140mm	19mm	44,5mm

Measures according to the DIN EN ISO 8624 standard.

QUEST CF

CALIBER



50 (12-14 years) 50mm 125-140mm 16mm 39,5mm
52 (14-16 years) 52mm 125-140mm 16mm 40,5mm

Measures according to the DIN EN ISO 8624 standard.



MATTE GREY/RED/BLUE

50: NAO3360250 52: NAO3360252



MATTE BLACK/BLUE TURQUOISE

50: NAO3360150 52: NAO3360152



MATTE DARK RASPBERRY/MAGENTA

50: NAO3360350 52: NAO3360352



CRYSTAL SOFT BLUE/BLUE

50: NAO3360450 52: NAO3360452




CRYSTAL PINK/DARK PINK


50: NAO3360550 52: NAO3360552



NANO

CUSTOM
FIT





NANO
CUSTOM
FIT



CALIBER    

52 (14-16 years)	52mm	125-140mm	15mm	38mm
54 (+16 years)	54mm	125-140mm	17mm	40mm

Measures according to the DIN EN ISO 8624 standard.



MATTE BLACK/GREEN
52: **NAO3370152** 54: **NAO3370154**



CRYSTAL NAVY/ORANGE
52: **NAO3370252** 54: **NAO3370254**



MATTE BLUE/GREY/ORANGE
52: **NAO3370452** 54: **NAO3370454**



MATTE DARK GREY/BLUE
52: **NAO3370352** 54: **NAO3370354**



CRYSTAL PEARL GREY/NAVY
52: **NAO3370552** 54: **NAO3370554**



NANO
VISTA_{m2}

$$\begin{aligned}\theta &= \omega_i(t_k - t_i) + \frac{1}{2} \frac{a(t_k - t_i)^2}{g \cos y} \\ a &= \frac{g \sin y}{1 + \frac{m}{M} \cos^2 y} \\ T &= 2\pi \sqrt{\frac{m}{k}} \quad E = k \frac{|Q|}{r^2} \\ F_s &= 6\pi\eta r v \\ \tau &= NIAB \sin \theta \\ \tan \theta &= \frac{v_z}{v_y}\end{aligned}$$

MYOP-K

STATE-OF-THE-ART
PROGRAM TO MANAGE
THE PROGRESSION OF
CHILDHOOD MYOPIA

FOR A
BRIGHTER LIFE!

“Each increasing diopter can be a lost opportunity”



Myopia can severely affect children's quality of life, hindering their learning, sports performance, and affecting their confidence. In the long term, it can lead to serious problems¹ such as:

- Myopic maculopathy
- Retinal detachment
- Cataracts
- Glaucoma

Delaying its progression by 1 diopter during childhood reduces by 40% the risk of myopic maculopathy, the leading cause of irreversible vision loss in the world².

1. Flitcroft DI. The complex interactions of retinal, optical and environmental factors in myopia etiology. *Prog Retin Eye Res.* 2012;31:622-60.
2. Bullimore MA, Brennan NA. Myopia Control: Why Each Diopter Matters. *Optom Vis Sci.* 2019 Jun;96(6):463-465

CHECK AVAILABILITY IN YOUR COUNTRY

THE PERFECT SYNERGY

The success of **MYOP-K** lies in its perfect synergy with **Nano Vista**, frames, forming an unbeatable combination. Our patented **Siliflex™** material is ultra-resistant, flexible, hypoallergenic, BPA-free, plant-based, and 100% recyclable, making our frames **#INDESTRUCTIBLE**.

Specifically designed for children, not just scaled-down adult glasses; we use children's facial scans for precise fitting at each growth stage.


Additionally, each frame has an interchangeable band and mini-band. This **versatility** allows adaptation to children's active lifestyle and keeps the lenses in the correct position all day.



NANO VISTA MYOP-K

HOW DOES IT WORK?

Nano Vista MYOP-K was born after years of research and development. Its revolutionary technology is based on a clear central zone and peripheral defocus that ensures all light rays focus within the retina, preventing greater ocular elongation and thus slowing myopia progression.



3. In 73% of children, myopia increase was less than 0.25D after one year of use.
4. Prospective, controlled, randomized, double-blind clinical study led by the European University of Madrid with 83 children (clinical trial NCT05250206).

EFFICACY DEMONSTRATED IN EUROPEAN STUDIES⁴

The results of the first study on European myopic children speak for themselves:

- In 73% of cases, MYOP-K lenses managed to significantly slow myopia progression compared to standard monofocal lenses.
- 91% of children using MYOP-K had an increase of less than 0.5D in both eyes over one year.

Additionally, MYOP-K stood out for its **comfort, performance, and clarity**, ensuring effective myopia treatment and excellent adaptability.

73% Myopia significantly slowed after 12 months



See complete study

NANO
VISTA^{my}

MYOP-K



Undetectable
aesthetics

WHY CHOOSE NANO VISTA MYOP-K?

- 1** Proven results in managing myopia progression.
- 2** Improves child's quality of life, both short and long term.
- 3** Excellent visual quality and quick adaptation period.
- 4** Aesthetically undetectable.
- 5** Perfect synergy with our Nano Indestructible frames.



Watch video

* The information in this document is general and does not constitute medical advice.
Consult a healthcare professional before using Nano Vista MYOP-K lenses.



nano CUSTOM
FIT



nano
VISTA
nano-vista.com

Phone: +34 600 514 183
email: info@gvo-optic.com

PRINT IN SPAIN

Gvo optics
gvo-optic.com

NANO
CUSTOM FIT
VIDEO



Authorized distributor:

CATN12EN

